

FOR ANIMAL TREATMENT ONLY

Porcilis APPvac®



Pleuropneumonia Vaccine for Pigs

ACTIVE CONSTITUENTS:

Inactivated adjuvanted vaccine containing:

Whole cell *Actinobacillus pleuropneumoniae* serotypes 1, 7 and 15 including a minimum 18 haemolysin units/mL

Preservatives: 0.1 mg/mL Thiomersal; < 0.05 mg/mL Formaldehyde

An aid in the control of pleuropneumoniae in pigs, caused by *Actinobacillus pleuropneumoniae* serotypes 1, 7 and 15.

DESCRIPTION:

Porcilis APPvac is a killed, whole cell vaccine containing *Actinobacillus pleuropneumoniae* (serovars 1, 7 and 15). It is recommended as an aid in control of pleuropneumonia in pigs. The adjuvant system increases the level of specific antibodies which act as an aid in protecting pigs through to market age.

SAFETY:

Porcilis APPvac has been shown to be safe when administered to pigs around 8 weeks of age at the recommended dose rate.

DIRECTIONS FOR USE:

DISCARD IF PREVIOUSLY FROZEN.

SHAKE WELL TO MIX BEFORE USE AND KEEP MIXED DURING USE.

Once opened, vaccine should be used within 24 hours.

Vaccinate pigs with 2mL by intramuscular injection at around 8 weeks of age and again 2-3 weeks later.

The recommended site of injection is that used for routine injection in pigs i.e. behind and below the ear.

CAUTION: Avoid Carcass Damage

1. Sterilise all injection apparatus by boiling before use. Avoid use of strong disinfectants on apparatus.
2. Maintain cleanliness at all times.
3. Keep needles sharp and clean. Replace frequently.
4. Use needles of appropriate gauge and length.
5. As far as possible avoid injection of animals during wet weather or under dusty conditions.
6. This product must be injected only into muscle tissue.
7. If possible inject into muscle tissue on side of neck.

CONNECTING PACK TO AUTOMATIC VACCINATOR:

Connect the tube firmly into the automatic vaccinator.

The pack may be suspended by means of a carrying strap.

Eject air from the vaccinator and tube by depressing plunger until vaccine flows through needle.

RESEALING VACCINE PACK:

For best practice vaccine should be used during one administration.

If however, unused vaccine is to be used within the following 24 hours, close clip as close as possible to the open end of the tube. Remove tube from the vaccinator gun.

Store between 2°C and 8°C. (Refrigerate. Do not freeze).

When re-using, cut tube between clip and vaccine pack and reconnect to the vaccinator gun.

Discard if previously frozen.

WITHHOLDING PERIOD: NIL

TRADE ADVICE:

EXPORT SLAUGHTER INTERVAL (ESI): ESI not required.

USER SAFETY AND FIRST AID INFORMATION:

The person who injects the vaccine should be careful to avoid accidental self-administration. Accidental self-administration may result in local bruising, pain and swelling. In the event of self-administration, seek medical attention and show the package leaflet or the label, to the Medical Practitioner.

If poisoning occurs, contact a Doctor or Poisons Information Centre. Phone Australia 131 126.

If swallowed do NOT induce vomiting. Give a glass of water.

DISPOSAL:

Dispose of empty container by wrapping with paper and putting in garbage.

STORAGE:

Store between 2°C and 8°C.

(Refrigerate. Do not freeze).

Protect from light.

Keep out of reach of children.



APVMA Approval No.: 63910/46889

MSD Animal Health
Intervet Australia Pty Limited
91-105 Harpin Street
BENDIGO EAST VIC 3550
Phone: 1800 033 461



MSD
Animal Health

141908


Client: MSD Animal Health

Product Description: Porcilis APPvac common Leaflet

Scale: 100%

Colour:  60% Black

 100% Black

Keyline:  (not to be printed)

Dimensions: 80 mm x 150 mm

Clients Code: 141908

PPC: 13-045

Agency Code & Version: MSD14891 V3

Last revised: 15/1/14

Signal Heading: Lucida Sans AK Bold 3 mm

Product Name: Lucida Sans AK Bold 4.6 mm

Body Text: Lucida Sans AK Roman 1.6 mm

Directions for Use: Lucida Sans AK Bold 1.6 mm

Safety Directions: Lucida Sans AK Bold 1.6 mm